MT CHAMBER SEAL \rightarrow LOCKING PROCEDURE




Diagram A




Diagram B




Diagram C




Diagram D

STEP 1

Detach each seal from the mat of 5 by first bending it backwards at the join followed by forward bending till it breaks, as shown in Diagram A

STEP 2

Hold the seal using one hand with the locking facing sideways, as shown in Diagram B

STEP 3

Insert the seal at the entry hole of the chamber, as shown in Diagram C.

STEP 4

After inserting, push the seal till you hear a click to indicate that it has securely locked, as shown in Diagram D.